


RATIO APPLICATOR KIT



The simple solution for simultaneous delivery of two liquids in a proven, effective ratio.

✓ **SIMPLE ONE-HANDED PRECISION**

Easily apply a uniform 12:1 mixture from two syringes with the single holder/plunger

✓ **FINELY MICRONIZED SPRAY**

Form a highly concentrated mist with an exceptionally small droplet size (20–230 microns)

✓ **ELITE PERFORMANCE STANDARDS**

Dispense 13mL in total deliverable volume with optimized flow rate, droplet size and velocity, and coverage

✓ **STEADY, ACCURATE STREAM**

Engineered to ensure uniform coverage through a consistent conical plume at the correct psi, optimized for sterile fluid application

Indication

The components of the Recibio, Inc. CeaLogic Ratio Applicator Kit are intended to assist the user in simultaneously delivering two liquids to the treatment site.

www.cealogic.com

The CeaLogic Ratio Applicator Kit is purpose-built for the simultaneous delivery of two fluids to the treatment site with simple, single-handed precision.

PRODUCT HIGHLIGHTS



A PROVEN RATIO

Designed for easy mixing from two, provided syringes at 12:1 ratio.




ANALYSIS-READY

The specimen cup and cap allow for sterile field transfer.



UNOBSTRUCTED VIEW

Transparent syringes offer a clear view and precise, calibrated measurement of fluids.



All CeaLogic components are latex free and are shipped sterile in package.

Qty.	Component
1	Sprayer Tip
1	Sprayer Dispenser
1	Sprayer Adapter
1	Specimen Cup
1	Specimen Label
1	Plastic Cap with Foam Liner
1	12cc Luer Lock Syringe
1	1cc Luer Lock Syringe
1	Product Sticker Sheet

