

SPECIMEN COLLECTION KIT



The all-in-one package for easy, safe collection and seamless, sterile transfer of fluid specimens.

✓ **SEAMLESS, STERILE COLLECTION**

Preserve sterility during transfer by drawing the specimen directly from the collection container

✓ **UNOBSTRUCTED SPECIMEN EVALUATION**

Clearly visualize and assess the specimen with the transparent collection container

✓ **EASY, INTUITIVE WORKFLOW**

Easily draw fluid from the specimen container's patented 45°-angled Luer lock™ port

✓ **INTEGRATED VOLUME MEASUREMENT**

Ensure precise measurement and maximum specimen collection with the 80cc calibrated container

Indication

The Recipro, Inc. CeaLogic Specimen Collection Kit is intended for normal suction and specimen collection.

www.cealogic.com

The CeaLogic Specimen Collection Kit's patented design makes it safe and simple to collect and transfer sterile fluid.

PRODUCT HIGHLIGHTS



SAFETY FIRST

Narrow, round-tip Yankauer suction promotes safe specimen collection.



ANALYSIS-READY

The specimen cup and cap allow for sterile field transfer.



FLEXIBLY FITS YOUR WORKFLOW

Extended tubing and 5-in-1 connector allow for connection to any vacuum.



All CeaLogic components are latex free and are shipped sterile in package.

Qty.	Component
1	Specimen Collection Container
1	10'-1/4" Suction Tubing with Male Adapter
1	Yankauer Suction Tip
1	Specimen Cap
1	Specimen Label
1	Product Sticker Sheet
1	5-in-1 Connector

